

AD-A155 736

LEPL

AD

AD-E401 352

TECHNICAL REPORT ARPAD-TR-85001

**AUTOMATIC SPHERICITY INTERFEROMETER
FOR TESTING LENS RADII**

JOHN SALERNO

LIBRARY

JUNE 1985



**U.S. ARMY ARMAMENT, MUNITIONS AND CHEMICAL COMMAND
PRODUCT ASSURANCE DIRECTORATE
DOVER, NEW JERSEY**

APPROVED FOR PUBLIC RELEASE; DISTRIBUTION UNLIMITED.

The views, opinions, and/or findings contained in this report are those of the author(s) and should not be construed as an official Department of the Army position, policy, or decision, unless so designated by other documentation.

The citation in this report of the names of commercial firms or commercially available products or services does not constitute official endorsement by or approval of the U.S. Government.

Destroy this report when no longer needed. Do not return to the originator.

REPORT DOCUMENTATION PAGE		READ INSTRUCTIONS BEFORE COMPLETING FORM
1. REPORT NUMBER Technical Report ARPAD-TR-85001	2. GOVT ACCESSION NO.	3. RECIPIENT'S CATALOG NUMBER
4. TITLE (and Subtitle) AUTOMATIC SPHERICITY INTERFEROMETER FOR TESTING LENS RADII		5. TYPE OF REPORT & PERIOD COVERED
		6. PERFORMING ORG. REPORT NUMBER
7. AUTHOR(s) John Salerno		8. CONTRACT OR GRANT NUMBER(s)
9. PERFORMING ORGANIZATION NAME AND ADDRESS ARDC, PAD Fire Control & Small Caliber Armament Systems Div (AMSMC-QAF-I) Dover, NJ 07801-5001		10. PROGRAM ELEMENT, PROJECT, TASK AREA & WORK UNIT NUMBERS QA Project M 6350 AMS 5397 OM 6350 AMCCOM MTT 50
11. CONTROLLING OFFICE NAME AND ADDRESS ARDC, TSD STINFO Div (SMCAR-TSS) Dover, NJ 07801-5001		12. REPORT DATE June 1985
		13. NUMBER OF PAGES 24
14. MONITORING AGENCY NAME & ADDRESS (If different from Controlling Office) Director U.S. Army Materials & Mechanics Research Center ATTN: DRXMR-STQ Watertown, MA 02172		15. SECURITY CLASS. (of this report) Unclassified
		15a. DECLASSIFICATION/DOWNGRADING SCHEDULE
16. DISTRIBUTION STATEMENT (of this Report) Approved for public release; distribution unlimited.		
17. DISTRIBUTION STATEMENT (of the abstract entered in Block 20, if different from Report)		
18. SUPPLEMENTARY NOTES This project has been accomplished as part of the U.S. Army's Manufacturing Testing Technology Program, which has for its objective the timely establishment of testing techniques, procedures, and prototype equipment (in mechanical, chemical, and nondestructive testing) to ensure efficient inspection (cont)		
19. KEY WORDS (Continue on reverse side if necessary and identify by block number) Interferometer Fringes Convex Digital radius scale Master radius standards Spherical wave Radii Test plates MTT-Radius measurements Transmission sphere Concave		
20. ABSTRACT (Continue on reverse side if necessary and identify by block number) The majority of lens radii measurements are made with test plates. This technique has the disadvantage of bringing the lens into direct contact with the test plate which can cause scratching of the surface and subsequent rejection of the test unit. In addition, test plate measurements are subjective because they depend on the experience and skill of the operator. As an alternative, an automated interferometric technique was evaluated for its ability to measure lens radii. This method uses a fizeau (cont)		

18. SUPPLEMENTARY NOTES (cont)

methods for material/materiel procured and maintained by AMC.

20. ABSTRACT (cont)

interferometer with the following accessories: transmission sphere, lens mount, and digital radius slide. This technique, chosen for its ease of operation and noncontact approach, proved to be a more suitable approach to radius measurement than the use of test plates, because it eliminated both the subjectivity and contact associated with the test plates.

CONTENTS

	Page
Introduction	1
Discussion	1
Conclusions	9
Recommendations	9
Distribution List	13

TABLES

1 Master test lenses	3
2 Calibrated steel balls	3
3 Interferometer readings of radius test standards	4
4 Interferometric readings of calibrated steel balls	8

FIGURES

1 Interferometer setup for measurement of concave lens radii (position 1)	11
2 Interferometer setup for measurement of concave lens radii (position 2)	11
3 Interferometer setup for measurement of convex lens radii (positions 1 and 2)	12
4 A graph of attainable system accuracy as a function of R/number	12

INTRODUCTION

This project involved the development and evaluation of an automated interferometer system for measuring the radius of curvature and surface figure of a spherical lens surface. An alternative to the use of conventional test glasses, the new interferometer system has the advantage of being a noncontacting technique. The conventional test glass technique involves physical contact between the test glass and the lens surface, a significant cause of lens scratching (particularly with inexperienced optical inspectors).

The equipment selected for this program was the Zygo MK11 Interferometer with a Digital Radius Scale (DRS) and a Zygo Automated Pattern Processor (ZAPP).¹ To evaluate the new system for accuracy, reliability, and ease of use, a program was devised to use the interferometer system to measure the curvature and figure of various lenses and compare the results with conventional measurements.

DISCUSSION

Before the interferometer is used, it must be set up with its optical axis aligned parallel to the DRS. Additionally, the accuracy of the measurements is based on the straightness of the guide bar and the flatness of the surface on which the 3-axis mount slides. In our case, the horizontal plane offset is 8 inches; therefore, an angular error of 1 arc second will introduce an error in the measured radius of 40 micro-inches. In the vertical plane the offset is 2 inches; therefore, an angular error of 1 arc second will introduce an error of 10 micro-inches. (Details of the alignment procedure are provided by the manufacturer of the interferometer system.^{2,3})

Samples of master radius test standards (table 1) obtained from an in-house optics shop were used to make radius measurements. The standards were selected to represent a cross section of lenses typically encountered in the Army's inventory. The optical shop foreman when asked about the technique used for the measurement of the standards replied that the masters do not need to be measured, which implied that the masters are perfect. Therefore the accuracy of the standards was not known, and there was no baseline to compare our measurements. For this reason, calibrated steel balls were obtained from the Metrology Laboratory (table 2).

¹ However, any other interferometer with the same associated equipment can be used to take radius measurements like those presented in this report.

² "Zygo Interferometer System Operation - Maintenance Manual and Warranty," Zygo Corporation, Middlefield, CT, February 1980.

³ "Zygo Radius Scale Operation - Maintenance Manual," Zygo Corporation, Middlefield, CT, April 1979.

Three radii of curvature measurements were made by different inspectors on each of the items in table 1, with the item under test being dismounted and re-mounted between measurements. In addition, three transmission spheres provided overlapping capability, and more than one was used to record the radius measurements. To select the optimum transmission sphere R/number to fill the aperture of a concave or convex surface, the R/number of the surface is calculated by the following formula:

$$R/\text{number} = \frac{\text{Radius of curvature of surface under test}}{\text{Clear aperture of surface under test}}$$

The transmission spheres used in this program were the f.75, f1.5 and the f3.3.

The procedure for making the radius measurements is straightforward and easy to learn without experience on the interferometer. An additional advantage of this method is that the test item does not come into contact with the transmission sphere. The procedures (used in this program) may have some elements peculiar to the Zygo interferometer, but in general they apply to any interferometer.

In the first step, the transmission sphere is placed in the accessory receptacle of the interferometer and aligned with the auto-align system. The test lens is mounted in the self-centering element holder so that the lens surface faces the interferometer aperture. The self-centering element holder is mounted in a 3-axis mount attached to the DRS carriage.

The radius of curvature of either a concave or convex surface can be measured. The center of curvature of the test surface is interferometrically made to coincide with the focus of the spherical wave emanating from the transmission sphere in figures 1 and 3 (position 1). The Z-axis fine adjustment screw is adjusted to optimize the straightness of the fringes. Final adjustment of this screw is to be in a clockwise direction and is assured by moving through best focus if necessary (position 1). The DRS display is cleared and the mount and carriage assembly are translated to place the test surface at the focus of the output beam from the interferometer transmission sphere in figures 2 and 3 (position 2). The Z-axis fine adjustment screw of the mount is adjusted again to optimize the straightness of the fringes. Final adjustment is to be in a clockwise direction and is assured by moving back through best focus if necessary (position 2). The radius of curvature of the surface is now displayed on the digital readout. The sign display will be plus for convex samples and minus for concave.

The accuracy to which the radius of curvature of a high quality spherical surface can be determined is a function of three things:

1. The R/number of the test surface
2. The accuracy of judging fringe straightness
3. The resolution and accuracy of the radius slide readout

The R/number is defined as the radius of curvature of the test surface divided by the clear aperture of the test surface.

The DRS display is a half-thousandth resolution system; each increment of measurement increases the display reading by 0.0005 inches.

Assuming the ability to judge straightness to within 1/10 fringe, the graph in figure 4 shows the accuracy of radius measurements as a function of the R/number of the surface under test. For example, an R1.5 surface can be measured to 1 μ m, and an R/7 surface to 25 μ m (0.0001 inch).

The results from the radius measurements on the master radius standards, listed by inspector and the R/number of the transmission sphere in table 3, show that the readings are repeatable and consistent between inspectors. The results also indicate a discrepancy of 34 thousandths of an inch in the worst case and between 2 and 4 ten-thousandths of an inch typically. Since the radius standards were an unknown, steel calibration balls were obtained to further assess the capability of the interferometer system to measure radii (table 4). The results, compared to the actual reading in table 2, confirm the accuracy and resolution of the interferometer system.

Table 1. Master test lenses

<u>Concave (in.)</u>	<u>Convex (in.)</u>
-0.574	+0.574
-1.050	+1.050
-1.651	+1.651
-3.232	+3.232
-4.275	+5.251
-5.251	+6.256
-6.256	
-14.262	
-19.330	

Table 2. Calibrated steel balls

<u>Nominal size (in.)</u>	<u>Measured size (in.)</u>
0.500	0.5001
0.625	0.62484
0.750	0.74997
0.875	0.87498
1.000	1.00003

Table 3. Interferometer readings of radius test standards

<u>Inspector</u>	<u>Transmission sphere</u>	<u>Test standard (in.)</u>	<u>Measured radii (in.)</u>
1	f.75	-3.232	-3.2355
			-3.2355
			-3.2310
		-4.275	-4.2760
			-4.2755
			-4.2755
		-5.251	-5.2550
			-5.2555
			-5.2555
	f1.5	-3.232	-3.2355
			-3.2355
			-3.2355
		+3.232	-3.2355
			-3.2350
			-3.2355
		-4.275	-4.2750
			-4.2750
			-4.2755
		-5.251	-5.2550
			-5.2550
			-5.2555
	f3.3	-3.232	-3.2355
			-3.2355
			-3.2350
		+3.232	+3.2360
			+3.2355
			+3.2355
		-4.275	-4.2750
			-4.2755
			-4.2750
		-5.251	-5.2545
			-5.2545
			-5.2555
		+5.251	+5.2525
			+5.2525
			+5.2525

Table 3. (cont)

<u>Inspector</u>	<u>Transmission sphere</u>	<u>Test standard (in.)</u>	<u>Measured radii (in.)</u>
2	f.75	-0.574	-0.5730
			-0.5735
			-0.5735
		+0.574	+0.5735
			+0.5735
			+0.5735
		-1.050	-1.0495
			-1.0495
			-1.0495
		+1.050	+1.0490
			+1.0490
			+1.0490
	f1.5	-1.651	-1.6520
			-1.6520
			-1.6520
		-1.050	-1.0490
			-1.0490
			-1.0495
		+1.050	-1.0490
			-1.0490
			-1.0495
		-0.574	-0.5735
			-0.5730
			-0.5730
		+0.574	+0.5730
			+0.5730
			+0.5730
		-1.651	-1.6525
			-1.6525
			-1.6525
		+1.651	+1.6520
			+1.6525
			+1.6525
		-6.256	-6.2530
			-6.2525
			-6.2530

Table 3. (cont)

<u>Inspector</u>	<u>Transmission sphere</u>	<u>Test standard (in.)</u>	<u>Measured radii (in.)</u>
		-14.262	-14.2295 -14.2300 -14.2290
		-19.330	-19.3240 -19.3245 -19.3245
	f3.3	-1.050	-1.0495 -1.0495 -1.0490
		+1.050	+1.0485 +1.0490 +1.0490
		-1.651	-1.6525 -1.6525 -1.6525
		+1.651	+1.6525 +1.6520 +1.6525
		-6.256	-6.2525 -6.2525 -6.2525
		+6.256	+6.2535 +6.2535 +6.2535
		-14.212	-14.2290 -14.2290 -14.2290
		-19.330	-19.3245 -19.3250 -19.3245
3	f.75	-3.232	-3.2360 -3.2360 -3.2355
		-4.275	-4.2760 -4.2755 -4.2755

Table 3. (cont)

<u>Inspector</u>	<u>Transmission sphere</u>	<u>Test standard (in.)</u>	<u>Measured radii (in.)</u>
		-5.251	-5.2555 -5.2555 -5.2560
	f1.5	-3.232	-3.2355 -3.2355 -3.2355
		+3.232	-3.2360 -3.2355 -3.2355
		-4.275	-4.2750 -4.2755 -4.2750
		-5.251	-5.2555 -5.2555 -5.2555
		-6.256	-6.2525 -6.2525 -6.2530
		-14.262	-14.2285 -14.2285 -14.2285
		-19.330	-19.3250 -19.3245 -19.3250
	f3.3	-3.232	-3.2355 -3.2360 -3.2360
		+3.232	+3.2350 +3.2355 +3.2355
		-6.256	-6.2525 -6.2525 -6.2525

Table 3. (cont)

<u>Inspector</u>	<u>Transmission sphere</u>	<u>Test standard (in.)</u>	<u>Measured radii (in.)</u>
		+6.256	+6.2540 +6.2540 +6.2540
		-14.262	-14.2285 -14.2275 -14.2280
		-19.330	-19.3250 -19.3250 -19.3250

Table 4. Interferometric readings of calibrated steel balls

<u>Transmission sphere</u>	<u>Nominal diameter (in.)</u>	<u>Measured radii (in.)</u>
f1.5	0.500	0.2500
		0.2500
		0.2500
	0.625	0.3125
		0.3125
		0.3125
	0.750	0.3750
		0.3750
		0.3750
	0.875	0.4375
		0.4375
		0.4375
	1.000	0.5000
		0.5000
		0.5000

CONCLUSIONS

This project proved that an interferometer and radius scale can be used to measure lens radii readily with little operator experience. The interferometric technique is a noncontact method of measurement which eliminates one of the major drawbacks of using test plates. An additional advantage of the interferometer is its reduced sensitivity to the level of skill of an operator, which avoids the problem associated with radius measurements using test plates. Test results show that the actual value of the test plates can introduce a sizable error.

The interferometer offers the additional benefit of being able to measure both the surface distortion and the radius of the test lens simultaneously. The interferometer can also be used to measure the surface figure of many other types of components (e.g., prisms, mirrors, and flats).

RECOMMENDATIONS

The interferometer is a worthwhile investment for any optical shop. Both versatile and easy to use, it eliminates the drawbacks associated with test plates, including the necessity to maintain a large inventory of various sizes. Test results show that the interferometer provides a more reliable measurement of radius than the test glasses. The only restrictions are the availability of transmission spheres and the length of the Digital Radius Scale used to make the measurements.




Figure 1. Interferometer setup for measurement of concave lens radii (position 1)




Figure 2. Interferometer setup for measurement of concave lens radii (position 2)




Figure 3. Interferometer setup for measurement of convex lens radii (positions 1 and 2)




Figure 4. A graph of attainable system accuracy as a function of R/number

DISTRIBUTION LIST

Metals and Ceramics Information Center
505 King Avenue
ATTN: Harold Mindlin, Director
James Lynch, Assistant Director
Columbus, OH 43201

Commander
Defense Technical Information Center
ATTN: Accessions Division (12)
Cameron Station
Alexandria, VA 22314

Commander
U.S. Army Foreign Science
and Technology Center
ATTN: DRXST-SD3
220 Seventh Street, N.E.
Charlottesville, VA 22901

Office of the Deputy Chief of Staff
for Research, Development, and Acquisition
ATTN: DAMA-ARZ-E
DAMA-CSS
Washington, DC 20310

Commander
U.S. Army Research Office
ATTN: George Mayer
J. J. Murray
P.O. Box 12211
Research Triangle Park, NC 27709

Commander
U.S. Army Materiel Development and
Readiness Command
ATTN: AMCQA-E
AMCQA-P
AMCDE-D
AMCDMD-FT
AMCLDC
AMCMT
AMCMM-M
5001 Eisenhower Avenue
Alexandria, VA 22333

Commander
U.S. Army Electronics Research
and Development Command
ATTN: DRSEL-PA-E, Stan Alster
J. Quinn
Fort Monmouth, NJ 07703

Commander
U.S. Army Missile Command
ATTN: DRSMI-TB (2)
DRSMI-TK, J. Alley
DRSMI-M
DRSMI-ET, Robert O. Black
DRSMI-QS, George L. Stewart, Jr.
DRSMI-EAT, R. Talley
DRSMI-QP
Redstone Arsenal, AL 35809

Commander
U.S. Army Materiel Systems
Analysis Activity
ATTN: DRXSY-MP
Aberdeen Proving Ground, MD 21005-5066

Director
Ballistic Research Laboratory
ATTN: AMXBR-OD-ST
Aberdeen Proving Ground, MD 21005-5066

Commander
U.S. Army Troop Support and
Aviation Materiel Readiness Command
ATTN: DRSTS-PL, J. Corwin (2)
DRSTS-Q
DRSTS-M
4300 Goodfellow Boulevard
St. Louis, MO 63120

Commander
U.S. Army Natick Research
and Development Command
ATTN: DRDNA-EM
Natick, MA 01760

Commander
U.S. Army Armament, Munitions and
Chemical Command
ATTN: AMSMC-QAF-I(D)
AMSMC-QAS(D), B. Aronowitz (5)
AMSMC-QAS-T(D), G. Zamloot (10)
Dover, NJ 07801-5001

Commander
U.S. Army Mobility Equipment Research
and Development Command

ATTN: DRDME-D
DRDME-E
DRDME-G
DRDME-H
DRDME-M
DRDME-T
DRDME-TQ
DRDME-V
DRDME-ZE
DRDME-N

Fort Belvoir, VA 22060

Commander
U.S. Army Tank-Automotive Materiel
Readiness Command

ATTN: DRSTA-Q (2)
Warren, MI 48090

Commander
U.S. Army Armament, Munitions
and Chemical Command

ATTN: AMSMC-QA(R) (2)
AMSMC-SC(R)
DRSAR-RDP(R)
AMSMC-EN(R)
AMSMC-QAE(R)
AMSMC-LEP-L(R)

Rock Island, IL 61299-6000

Commander
Rock Island Arsenal
ATTN: SMCRI-EN, W. M. Kisner
SMCRI-ENM, W. D. McHenry
SMCRI-QA

Rock Island, IL 61299

Commander
U.S. Armament Research and Development Center
U.S. Army Armament, Munitions
and Chemical Command

ATTN: SMCAR-LC, E. Kelly
SMCAR-LCA
SMCAR-LCE
SMCAR-SCM, J. D. Corrie
SMCAR-SCM-O, Harry E. Pebly, Jr.
SMCAR-TSP, B. Stephans
SMCAR-TSS (5)

Dover, NJ 07801

Commander
U.S. Army Armament, Munitions and
Chemical Command
ATTN: AMSMC-GCL(D)
Dover, NJ 07801-5001

Commander
Chemical Research and Development Center
U.S. Army Armament, Munitions
and Chemical Command
ATTN: SMCCR-SPS-I
Aberdeen Proving Ground, MD 21010-5423

Commander
Chemical Research and Development Center
U.S. Army Armament, Munitions
and Chemical Command
ATTN: SMCCR-RSP-A
Aberdeen Proving Ground, MD 21010-5423

Commander
Armament Research and Development Center
U.S. Army Armament, Munitions
and Chemical Command
Product Assurance Directorate
ATTN: SMCAR-QAC-E(A), W. J. Maurits
Aberdeen Proving Ground, MD 21010

Chief
Benet Weapons Laboratory, LCWSL
Armament Research and Development Center
U.S. Army Armament, Munitions
and Chemical Command
ATTN: SMCAR-LCB, T. Moraczewski
SMCAR-LCB-TL
SARWV-PPI, L. Jette
Watervliet, NY 12189-5000

Commander
U.S. Army Aviation Research
and Development Command
ATTN: DRDAV-EXT
DRDAV-QR
DRDAV-QP
DRDAV-QE
St. Louis, MO 63120

Commander
U.S. Army Tank-Automotive Research
and Development Command
ATTN: DRDTA-UL, Technical Library
DRDTA-RCKM, S. Goodman
DRDTA-RCKT, J. Fix
DRDTA-RTAS, S. Catalano
DRDTA-TTM, W. Moncrief
DRDTA-ZS, O. Renius
DRDTA-JA, C. Kedzior
Warren, MI 48090

Director
Industrial Base Engineering Activity
ATTN: DRXIB-MT
Rock Island, IL 61299

Commander
Harry Diamond Laboratories
ATTN: DELHD-EDE, B. F. Willis
2800 Powder Mill Road
Adlephi, MD 20783

Commander
U.S. Army Test and Evaluation Command
ATTN: DRSTE-TD
DRSTE-ME
Aberdeen Proving Ground, MD 21005

Commander
U.S. Army White Sands Missile Range
ATTN: STEWS-AD-L
STEWS-ID
STEWS-TD-PM
White Sands Missile Range, NM 88002

Commander
U.S. Army Yuma Proving Ground
ATTN: Technical Library
Yuma, AZ 85364

Commander
U.S. Army Tropic Test Center
ATTN: STETC-TD, Drawer 942
Fort Clayton, Canal Zone

Commander
Aberdeen Proving Ground
ATTN: STEAP-MT
STEAP-MT-M, J. A. Feroli
STEAP-MT-G, R. L. Huddleston
Aberdeen, MD 21005

Commander
U.S. Army Cold Region Test Center
ATTN: STECR-OP-PM
APO Seattle 98733

Commander
U.S. Army Dugway Proving Ground
ATTN: STEDP-MT
Dugway, UT 84022

Commander
U.S. Army Electronic Proving Ground
ATTN: STEEP-MT
Fort Huachuca, AZ 35613

Commander
Jefferson Proving Ground
ATTN: STEJP-TD-I
Madison, IN 47250

Commander
U.S. Army Aircraft Development
Test Activity
ATTN: STEBG-TD
Fort Rucker, AL 36362

President
U.S. Army Armor and Engineer Board
ATTN: ATZKOE-TA
Fort Knox, KY 40121

President
U.S. Army Field Artillery Board
ATTN: ATZR-BDOP
Fort Sill, OK 73503

Commander
Anniston Army Depot
ATTN: SDSAN-QA
Anniston, AL 36202

Commander
Corpus Christi Army Depot
ATTN: SDSCC-MEE, Mr. Haggerty, Mail Stop 55
Corpus Christi, TX 78419

Commander
Letterkenny Army Depot
ATTN: SDSLE-QA
Chambersburg, PA 17201

Commander
Lexington-Bluegrass Army Depot
ATTN: SDSLX-QA
Lexington, KY 40507

Commander
New Cumberland Army Depot
ATTN: SDSNC-QA
New Cumberland, PA 17070

Commander
U.S. Army Depot Activity
ATTN: SDSTE-PU-Q (2)
Pueblo, CO 81001

Commander
Red River Army Depot
ATTN: SDSRR-QA
Texarkana, TX 75501

Commander
Sacramento Army Depot
ATTN: SDSSA-QA
Sacramento, CA 95813

Commander
Savanna Army Depot Activity
ATTN: SDSSV-S
Savanna, IL 61074

Commander
Seneca Army Depot
ATTN: SDSSE-R
Romulus, NY 14541

Commander
Sharpe Army Depot
ATTN: SDSSH-QE
Lathrop, CA 95330

Commander
Sierra Army Depot
ATTN: SDSSI-DQA
Herlong, CA 96113

Commander
Tobyhanna Army Depot
ATTN: SDSTO-Q
Tobyhanna, PA 18466

Commander
Tooele Army Depot
ATTN: SDSTE-QA
Tooele, UT 84074

Director
DARCOM Ammunition Center
ATTN: SARAC-DE
Savanna, IL 61074

Naval Research Laboratory
ATTN: Code 5830, J. M. Krafft
Code 2620, Library
Washington, DC 20375

Air Force Materials Laboratory
Wright-Patterson Air Force Base
ATTN: AFML-LTM, W. Wheeler
AFML-LLP, R. Rowand
Wright-Patterson Air Force Base, OH 45433

Director
U.S. Army Materials and
Mechanics Research Center
ATTN: DRXMR-D
DRXMR-P
DRXMR-S
DRXMR-STQ (4)
DRXMR-PL (2)
DRXMR-X
DRXMR-PD
DRXMR-WD
DRXMR-M (2)
DRXMR-ST
DRXMR-L
DRXMR-T
DRXMR-E
DRXMR-PR

Watertown, MA 02172

Columbia Research Corporation
1 Metropolitan Grove Road
ATTN: David Flake (2)
George Fisher (2)
Paul Watson (2)
Gaithersburg, MD 20760